

## Product datasheet for **TP304984L**

### ARF4 (NM\_001660) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ADP-ribosylation factor 4 (ARF4), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204984 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MGLTISSLFSRLFGKKQMRILMVGLDAAGKTTILYKLLKLGIVTTIPTIGFNVETVEYKNICFTVWDVGG QDRIRPLWKHYFQNTQGLIFVDSNDRERIQEVADELQKMLLVDELRAVLLLFANKQDLPNAMAISEMT DKLGLQSLRNRTWYVQATCATQGTGLYEGLDWLSNELSKR  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	20.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001651</a>
Locus ID:	378
UniProt ID:	<a href="#">P18085</a>
RefSeq Size:	1758



[View online »](#)


Cytogenetics: 3p14.3

RefSeq ORF: 540

Synonyms: ARF2

**Summary:** This gene is a member of the human ARF gene family whose members encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking and as activators of phospholipase D. The gene products include 5 ARF proteins and 11 ARF-like proteins and constitute one family of the RAS superfamily. The ARF proteins are categorized as class I, class II and class III; this gene is a class II member. The members of each class share a common gene organization. The ARF4 gene spans approximately 12kb and contains six exons and five introns. This gene is the most divergent member of the human ARFs. Conflicting map positions at 3p14 or 3p21 have been reported for this gene. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified ARF4 protein (Cat# [TP304984]). The protein was produced from HEK293T cells transfected with ARF4 cDNA clone (Cat# [RC204984]) using MegaTran 2.0 (Cat# [TT210002]).