

Product datasheet for TP305005M

P2Y12 (P2RY12) (NM_022788) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human purinergic receptor P2Y, G-protein coupled, 12 (P2RY12), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205005 protein sequence Red=Cloning site Green=Tags(s)

MQAVDNLTSAPGNTSLCTRDKYKITQVLFPLLYTVLFFVGLITNGLAMRIFFQIRSKSNFIIFLKNTVISD
LLMILTFPFKILSDAKLGTGPLRFTVCQVTSVIFYFTMYISISFLGLITIDRYQKTTRPFKTSNPKNLLG
AKILSVVIWAFMFLLSLPNMILTNRQPRDKNVKKCSFLKSEFGLVWHEIVNYICQVIFWINFLIVICYT
LITKELYRSYVRTRGVGKVPKRVNVKVFIIIVFFICFVPPHFARIPYTLTSQTRDVFDCCTAENTLFFVK
ESTLWLTSLNACLDPFIYFFLCKSFRNSLISMLKCPNSATLSLQDNRRKKEQDGGDPNEETPM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	39.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_073625
Locus ID:	64805



[View online »](#)

UniProt ID:	Q9H244 , A8K7T1
RefSeq Size:	2318
Cytogenetics:	3q25.1
RefSeq ORF:	1026
Synonyms:	ADPG-R; BDPLT8; HORK3; P2T(AC); P2Y(12)R; P2Y(AC); P2Y(ADP); P2Y(cyc); P2Y12; SP1999
Summary:	The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor is involved in platelet aggregation, and is a potential target for the treatment of thromboembolisms and other clotting disorders. Mutations in this gene are implicated in bleeding disorder, platelet type 8 (BDPLT8). Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, Jul 2013]
Protein Families:	Druggable Genome, GPCR, Transmembrane

Product images:



Coomassie blue staining of purified P2RY12 protein (Cat# [TP305005]). The protein was produced from HEK293T cells transfected with P2RY12 cDNA clone (Cat# [RC205005]) using MegaTran 2.0 (Cat# [TT210002]).