

Product datasheet for **TP305023**

CCDC12 (NM_144716) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 12 (CCDC12), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205023 protein sequence Red =Cloning site Green =Tags(s)
	MEATTAGVGRLEEEALRRKERLKAALREKTGRKDKEDGEPKTKHLREEEEEGEKHRELRLRNYVPEDEDLK KRRVPQAKPVAVEEKVKEQLEAAKPEPVIEEVDLANLAPRKPDPDWDLKRDVAKKLEKLEKRTQRAIAELIR ERLKGQEDSLASAVDAATEQKTCDS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_653317
Locus ID:	151903
UniProt ID:	J3KR35
RefSeq Size:	1035



[View online »](#)

Cytogenetics: 3p21.31
RefSeq ORF: 498
Protein Pathways: Spliceosome

Product images:

Coomassie blue staining of purified CCDC12 protein (Cat# TP305023). The protein was produced from HEK293T cells transfected with CCDC12 cDNA clone (Cat# [RC205023]) using MegaTran 2.0 (Cat# [TT210002]).