

## Product datasheet for TP305051L

### HAO2 (NM\_016527) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hydroxyacid oxidase 2 (long chain) (HAO2), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205051 protein sequence Red=Cloning site Green=Tags(s)

MSLVCLTDFQAHAREQLSKSTRDFIEGGADDSITRDDNIAAFKRIRLRPRYL RDVSEVDTRTTIQGEEIS  
APICIAPTGFHCLVWPDGEMSTARAAQAAGICYITSTFASCSLEDIVIAAPEGLRWFQLYVHPDLQLNKQ  
LIQRVESLGFKALVITLDT PVCGNRRHDIRNQLRRNLTLTDLQSPKKGNAIPYFQMTPISTSLCWNDSLW  
FQSITRLPIILKGILTKEDAELAVKHNVQGIIVSNHGGRQLDEVLASIDALTEVVA AVKKGKIEVYLDGGV  
RTGNDV LKALALGAKCIFLGRPILWGLACKGEHGVKEVLNILTNEFHTSMALTGCRSVAEINRNLVQFSR  
L

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	38.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_057611</a></u>



[View online »](#)

Locus ID:	51179
UniProt ID:	<a href="#">Q9NYQ3</a>
RefSeq Size:	1488
Cytogenetics:	1p12
RefSeq ORF:	1053
Synonyms:	GIG16; HAOX2
Summary:	This gene is one of three related genes that have 2-hydroxyacid oxidase activity. The encoded protein localizes to the peroxisome has the highest activity toward the substrate 2-hydroxypalmitate. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]
Protein Pathways:	Glyoxylate and dicarboxylate metabolism, Metabolic pathways

### Product images:



Coomassie blue staining of purified HAO2 protein (Cat# [TP305051]). The protein was produced from HEK293T cells transfected with HAO2 cDNA clone (Cat# [RC205051]) using MegaTran 2.0 (Cat# [TT210002]).