

Product datasheet for **TP305057L**

FAM151A (NM_176782) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 151, member A (FAM151A), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205057 protein sequence Red =Cloning site Green =Tags(s)

MVCREQLSKNQVKWVFAGITCVSVVIAAIVLAITLRRPGCELEACSPDADMLDYLLSLGQISRRDALEV
TWYHAANSKKAMTAALDSNITVLEADVNVVEGLGTANETGVPIMAHPPTTYSNDNTLEQWLDVAVLGSSQKGI
KLDKFNIAVGPDLRLRQLTEEGKVRRIWINADILKGNMLISTEVNATQFLALVQEYKPKATLSPGW
TTFYMSTSPNRTYTQAMVEKMHELVEGVPQVRVTFPVRSSMVRAAWPHFSWLLSQSERYSLTLWQAASDPM
SVEDLLYVRDNTAVHQVYYDIFELLSQFKLALNATRKPYYTGGSLIPLLQLPGDDGLNVEWLVDPVQ
GSGKTATMTLPDTEGMILLNTGLEGTVAENPVPIVHTPSGNILTLESCLQQLATHPGHWGIHLQIVEPAA
LRPSLALLARLSSLGLLHWPVWVGAKISHGFSVPGHVAGRELLTAVAEVFPHTVAPGWPEEVLGSGYR
EQLLTDMLELCQGLWQPVSFQMAMLLGHSTAGAIGRLLASSPRATVTVEHNPAGGDYASVRTALLAARA
VDRTRVYYRLPQGYHKDLLAHVGRN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	63.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_788954
Locus ID:	338094
UniProt ID:	Q8WW52
RefSeq Size:	1998
RefSeq ORF:	1755
Synonyms:	C1orf179
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified FAM151A protein (Cat# [TP305057]). The protein was produced from HEK293T cells transfected with FAM151A cDNA clone (Cat# [RC205057]) using MegaTran 2.0 (Cat# [TT210002]).