

Product datasheet for **TP305174L**

COQ10B (NM_025147) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coenzyme Q10 homolog B (<i>S. cerevisiae</i>) (COQ10B), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205174 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAARTGHTALRRVVGCRPKSATAAGAQAQAVRNGRYLASCILMSRTLPLHTSILPKEICARTFFKITAP LINKRKEYSERRILGYSMQEMYDVVSGVEDYKHFPWCKKSDVISKRSGYCKTRLEIGFPPVLERYTSV TLVKPHLVKASCTDGRLFNHLETIWRFPGLPGYPRTCTLDFSISFEFRSLLHSQATLFFDEVVKQMA AFERRACKLYGPETNIPRELMLHEVHHT</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	26.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_079423</u>
Locus ID:	80219
UniProt ID:	<u>Q9H8M1</u>



[View online »](#)


RefSeq Size: 1989

Cytogenetics: 2q33.1

RefSeq ORF: 714

Summary: Required for the function of coenzyme Q in the respiratory chain. May serve as a chaperone or may be involved in the transport of Q6 from its site of synthesis to the catalytic sites of the respiratory complexes (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified COQ10B protein (Cat# [TP305174]). The protein was produced from HEK293T cells transfected with COQ10B cDNA clone (Cat# [RC205174]) using MegaTran 2.0 (Cat# [TT210002]).