

Product datasheet for TP305249L

ARTS1 (ERAP1) (NM_001040458) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human endoplasmic reticulum aminopeptidase 1 (ERAP1), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205249 protein sequence Red =Cloning site Green =Tags(s)
	<p>MVFLPLKWSLATMSFLLSLLALLTVSTPSWCQSTEASPKRSDGTPFPWPNKIRLPEYVIPVHYDLLIHAN LTTLTFWGTTKVEITASQPTSTIILHSHHLQISRATLRKGAGERLSEEPLQVLEHPRQEQIALLAPEPLL VGLPYTVWIHYAGNLSETFHGFYKSTYRTKEGELRILASTQFEPTAARMAFPCFDEPAFKASFSIKIRRE PRHLAISNMPLVKSVTVAEGLIEDHFDVTVKMSTYLVAFIISDFESVSKITKSGVKVSVYAVPDKINQAD YALDAAVTLLEFYEDYFSIPYPLPKQDLAIPDFQSGAMENWGLTTYRESALLFDAEKSSASSKLGITMT VAHELAHQWFGNLVTMEWWNDLWLNNEGFAKFMFVSVSVTHPELKVGDYFFGKCFDAMEVDALNSSHPVS TPVENPAQIREMFDDVSYDKGACILNMLREYLSADAFKSGIVQYLQKHSYKNTKNEDLWDSMASICPTDG VKGMDGFCRSRQHSSSSSHHWHQEGVDVKTMMNTWTLQKGFPLITIVRGRNVHMKQEHYMKGSDGAPDTG YLWHVPLTFITSKSDMVHRFLLKTKTDVLILPEEVEWIKFNVGMNGYIVHYEDDGWDSLTGLLKGTHTA VSSNDRASLINNAFQLVLSIGKLSIEKALDLSLYLKHETEIMPVFQGLNELIPMYKLMKRDMEVETQFK AFLIRLLRDLIDKQWTWDEGSVSRMLRSQLLLLACVHNYQPCVQRAEGYFRKWKESNGNLSLPVDVTLA VFAVGAQSTEGWDFLYSKYQFSLSSSTEKSQIEFALCRTQNKEKLQWLLDESFKGDKIKTQEFQILTIG RNPVGYPLAWQFLRKNWNKLVQKFELGSSSIAHMVMGTTNQFSTRTRLEEVKGFSSLKENGSQLRCVQQ TIETIEENIGWMDKNFDKIRVWLQSEKLERM</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	107.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001035548
Locus ID:	51752
UniProt ID:	Q9NZ08
RefSeq Size:	5085
Cytogenetics:	5q15
RefSeq ORF:	2823
Synonyms:	A-LAP; ALAP; APPILS; ARTS-1; ARTS1; ERAAP; ERAAP1; PILS-AP; PILSAP
Summary:	The protein encoded by this gene is an aminopeptidase involved in trimming HLA class I-binding precursors so that they can be presented on MHC class I molecules. The encoded protein acts as a monomer or as a heterodimer with ERAP2. This protein may also be involved in blood pressure regulation by inactivation of angiotensin II. Three transcript variants encoding two different isoforms have been found for this gene.[provided by RefSeq, Oct 2010]
Protein Families:	Druggable Genome, Protease, Secreted Protein

Product images:



Coomassie blue staining of purified ERAP1 protein (Cat# [TP305249]). The protein was produced from HEK293T cells transfected with ERAP1 cDNA clone (Cat# [RC205249]) using MegaTran 2.0 (Cat# [TT210002]).