

Product datasheet for TP305271

RASL10A (NM_006477) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RAS-like, family 10, member A (RASL10A), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205271 protein sequence Red=Cloning site Green=Tags(s)

MGGSLRVAVLGAPGVGKTAIIRQFVFGDYPERHRPTDGPRLYRPAVLLDGAVYDLSIRDGDVAGPGSSPG
GPEEWPDAKDWSLQDTDAFVLVYDICSPDSFDYVKALRQRIAETRPAGAPEAPILVVGNGKRDRQRLRFGP
RRALAALVRRGWRCGYLECSAKYNWHVLRFLRELLRCALVRARPAHPALRLQGALHPARCSLM


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	22.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006468
Locus ID:	10633
UniProt ID:	Q92737 , A0A024R1C8



[View online »](#)

RefSeq Size:	1504
Cytogenetics:	22q12.2
RefSeq ORF:	609
Synonyms:	RRP22
Summary:	Potent inhibitor of cellular proliferation.[UniProtKB/Swiss-Prot Function]

Product images:

Coomassie blue staining of purified RASL10A protein (Cat# TP305271). The protein was produced from HEK293T cells transfected with RASL10A cDNA clone (Cat# [RC205271]) using MegaTran 2.0 (Cat# [TT210002]).