

Product datasheet for **TP305299**

GNG4 (NM_004485) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human guanine nucleotide binding protein (G protein), gamma 4 (GNG4), transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205299 protein sequence Red =Cloning site Green =Tags(s) MKEGMSNNSTTSISQARKAVEQLKMEACMDRVKVSQAAADLLAYCEAHVREDPLIIPVPASENPFREKKF FCTIL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	8.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004476
Locus ID:	2786
UniProt ID:	P50150 , B1APZ0
RefSeq Size:	4885



[View online »](#)

Cytogenetics: 1q42.3


RefSeq ORF: 225

Summary: Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction. [UniProtKB/Swiss-Prot Function]

Protein Families: Druggable Genome

Protein Pathways: Chemokine signaling pathway

Product images:



Coomassie blue staining of purified GNG4 protein (Cat# TP305299). The protein was produced from HEK293T cells transfected with GNG4 cDNA clone (Cat# [RC205299]) using MegaTran 2.0 (Cat# [TT210002]).