

## Product datasheet for TP305311L

### TULP2 (NM\_003323) Human Recombinant Protein

#### Product data:

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant protein of human tubby like protein 2 (TULP2), 1 mg
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>RC205311 protein sequence Red=Cloning site Green=Tags(s)
	MSQDNDTLMRDILGHELAAMRLQKLEQQRRLFEKKQRQKRQELLMVQANPDASPWLWRSCLREERLLGDR GLGNPFLRKKVSEAHLPSTGIHLSALGTVSCGGDGRGERGLPTPRTEAVFRNLGLQSPFLSWLPDNSDAELE EVSVENGSVPPPFKQSPRIRKKGWQAHQRPQTRAEGESDSQDMGDAHKSPNMGPNPGMDGDCVYENLAF QKEEDLEKKREASESTGTNSSAAHNEELSKALKGEGGTSDSDHMRHEASLAIRSPCPGLEEDMEAYVLRPA LPGTMMQCYLTRDKHGVKGLFPLYLYLETSDSLQRFLLAGRKRKRRSNTSYLISLDPHTLSRDGDNFV GKVRSNVFSTKFTIFDNGVNPDRHLTRNTARIRQELGAVCYEPNVLGYLGPRKMTVILPGTNSQNQRIN VQPLNEQESLLSRYQRGDKQGLLLLHNKTPSWDKENGVYTLNFHGRVTRASVKNFQIVDPKHQEHVLVLF GRVGPDTFTMDFCFPFSPQLQAFSICLSSFN
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
<b>Tag:</b>	C-Myc/DDK
<b>Predicted MW:</b>	58.5 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_003314](#)

Locus ID: 7288

UniProt ID: [O00295](#)

RefSeq Size: 1791


Cytogenetics: 19q13.33

RefSeq ORF: 1560

Synonyms: CT65; TUBL2

**Summary:** TULP2 is a member of a family of tubby-like genes (TULPs) that encode proteins of unknown function. Members of this family have been identified in plants, vertebrates, and invertebrates. The TULP proteins share a conserved C-terminal region of approximately 200 amino acid residues. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified TULP2 protein (Cat# [TP305311]). The protein was produced from HEK293T cells transfected with TULP2 cDNA clone (Cat# [RC205311]) using MegaTran 2.0 (Cat# [TT210002]).