

Product datasheet for **TP305364**

EBAG9 (NM_004215) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human estrogen receptor binding site associated, antigen, 9 (EBAG9), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205364 protein sequence Red =Cloning site Green =Tags(s) MAITQFRLFKFCTCLATVFSFLKRLICRSGRGRKLSGDQITLPTTVDYSSVPKQTDVEEWTSWDEDAPTS VKIEGNGNVATQQNSLEQLPFDYFKDMTPTIRKTQKIVIKKREPLNFGIPDGSTGFSSRLAATQDLPI HQSSGLDLDTWQENTNAWEEEEEDAAWQAEVLRQQKLADREKRAAEQQRKKMEKEAQRMLMKKEQNKIGV KLS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004206
Locus ID:	9166



[View online »](#)

UniProt ID: [O00559](#), [A0A024R9E0](#)

RefSeq Size: 2495

Cytogenetics: 8q23.2

RefSeq ORF: 639

Synonyms: EB9; PDAF

Summary: This gene was identified as an estrogen-responsive gene. Regulation of transcription by estrogen is mediated by estrogen receptor, which binds to the estrogen-responsive element found in the 5'-flanking region of this gene. The encoded protein is a tumor-associated antigen that is expressed at high frequency in a variety of cancers. Alternate splicing results in multiple transcript variants. A pseudogene of this gene has been defined on chromosome 10. [provided by RefSeq, Jul 2013]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified EBAG9 protein (Cat# TP305364). The protein was produced from HEK293T cells transfected with EBAG9 cDNA clone (Cat# [RC205364]) using MegaTran 2.0 (Cat# [TT210002]).