

Product datasheet for **TP305407**

RNF113B (NM_178861) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human ring finger protein 113B (RNF113B), 20 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC205407 protein sequence
Red=Cloning site **Green**=Tags(s)

MAAPPSPGRTADQADQVCTFLFKKPKGRKGAAGLRKRPACDPEHGESSSSGDEGDTVAQPPRVAPRPRGLH
SWQKAAHGDRRGEEAAPESLDWVYRSTRSAKPVGPEDMGATADFEQDTEKEHHTPTILKCSQRVQEALRG
REHDHIYRGIHSYLRYLKPKDTSMGNSSSGMARKGPIRAPGHLRATVRWDYQPDICKDYKETGFCGFGDS
CKFLHDRSDYKLGWEIERELEEGRYCICEDENHEVGSEEEEIPFRFCICRQAFQNPVVTCKRHYFCESCA
LEHFRATPRCYICDQPTGGIFNPAKELMAKLQKLQAAEGKKR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 36.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_849192](#)


Locus ID: 140432



[View online »](#)

UniProt ID:	Q8IZP6
RefSeq Size:	1381
Cytogenetics:	13q32.2
RefSeq ORF:	966
Synonyms:	bA10G5.1; RNF161; ZNF183L1
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified RNF113B protein (Cat# TP305407). The protein was produced from HEK293T cells transfected with RNF113B cDNA clone (Cat# [RC205407]) using MegaTran 2.0 (Cat# [TT210002]).