

Product datasheet for **TP305410L**

CAPSL (NM_001042625) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human calcyphosine-like (CAPSL), transcript variant 2, 1 mg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC205410 protein sequence
Red=Cloning site **Green**=Tags(s)

MAIQAKKKLTTATNPIERLRLQCLARGSAGIKGLGRVFRIMDDDDNNRTLDFKEFMKGLNDYAVVMEKEEV
EELFRRFDKDGNGTIDFNEFLTLRPPMSRARKEVIMQAFRKLDKTGDGVITIEDLREVYNAKHHPKYQN
GEWSEEQVFRKFLDNFDSPYDKDGLVTPPEEFMNYAGVSASIDTDVYFIIMMRTAWKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 24 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001036090](#)

Locus ID: 133690

UniProt ID: [Q8WWF8](#)

RefSeq Size: 987




[View online »](#)

Cytogenetics: 5p13.2

RefSeq ORF: 594

Product images:



Coomassie blue staining of purified CAPSL protein (Cat# [TP305410]). The protein was produced from HEK293T cells transfected with CAPSL cDNA clone (Cat# [RC205410]) using MegaTran 2.0 (Cat# [TT210002]).