

Product datasheet for **TP305453M**

HORMAD2 (NM_152510) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human HORMA domain containing 2 (HORMAD2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205453 protein sequence Red =Cloning site Green =Tags(s)
	<p>MATAQLSHCITIHKASKETVFPSPQITNEHESLKMVKKLFATSISCITYLRGLFPESYGERHLLDLSLKI LREDDKKCPGSLHIIRWIQGCFDALEKRYLRMAVLTLYTDPMGSEKVTEMYQKFKYTKEGATMDFDSHSS STSFESGTNNEDIKKASVLLIRKLYILMQDLEPLPNNVLTMKLHYNAVTPHDYQPLGFKGVNSHFL FDKEPINVQGVFSTGFHSMKVVMTEATKVIDLENNLFRENSTTEIAHQGLDCDEEEECNDHIQRMNFV CSQQSSECSRKKRKVSEPVKVFIPNRK</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	35.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_689723</u>
Locus ID:	150280


[View online »](#)

UniProt ID: Q8N7B1

RefSeq Size: 1909


Cytogenetics: 22q12.2

RefSeq ORF: 921

Synonyms: CT46.2

Summary: Essential for synapsis surveillance during meiotic prophase via the recruitment of ATR activity. Plays a key role in the male mid-pachytene checkpoint and the female meiotic prophase checkpoint: required for efficient build-up of ATR activity on unsynapsed chromosome regions, a process believed to form the basis of meiotic silencing of unsynapsed chromatin (MSUC) and meiotic prophase quality control in both sexes. Required for the DNA double-strand break-independent, BRCA1-dependent activation of ATR on the sex chromosomes that is essential for normal sex body formation (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified HORMAD2 protein (Cat# [TP305453]). The protein was produced from HEK293T cells transfected with HORMAD2 cDNA clone (Cat# [RC205453]) using MegaTran 2.0 (Cat# [TT210002]).