

Product datasheet for TP305473M

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

RASA1 (NM_002890) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human RAS p21 protein activator (GTPase activating protein) 1

(RASA1), transcript variant 1, 100 μg

Species: Human Expression Host: HEK293T

Expression cDNA Clone >RC205473 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MMAAEAGSEEGGPVTAGAGGGGAAAGSSAYPAVCRVKIPAALPVAAAPYPGLVETGVAGTLGGGAALGSE FLGAGSVAGALGGAGLTGGGTAAGVAGAAAGVAGAAVAGPSGDMALTKLPTSLLAETLGPGGGFPPLPPP PYLPPLGAGLGTVDEGDSLDGPEYEEEEVAIPLTAPPTNQWYHGKLDRTIAEERLRQAGKSGSYLIRESD RRPGSFVLSFLSQMNVVNHFRIIAMCGDYYIGGRRFSSLSDLIGYYSHVSCLLKGEKLLYPVAPPEPVED RRRVRAILPYTKVPDTDEISFLKGDMFIVHNELEDGWMWVTNLRTDEQGLIVEDLVEEVGREEDPHEGKI WFHGKISKQEAYNLLMTVGQVCSFLVRPSDNTPGDYSLYFRTNENIQRFKICPTPNNQFMMGGRYYNSIG DIIDHYRKEQIVEGYYLKEPVPMQDQEQVLNDTVDGKEIYNTIRRKTKDAFYKNIVKKGYLLKKGKGKRW KNLYFILEGSDAQLIYFESEKRATKPKGLIDLSVCSVYVVHDSLFGRPNCFQIVVQHFSEEHYIFYFAGE TPEQAEDWMKGLQAFCNLRKSSPGTSNKRLRQVSSLVLHIEEAHKLPVKHFTNPYCNIYLNSVQVAKTHA REGQNPVWSEEFVFDDLPPDINRFEITLSNKTKKSKDPDILFMRCQLSRLQKGHATDEWFLLSSHIPLKG IEPGSLRVRARYSMEKIMPEEEYSEFKELILQKELHVVYALSHVCGQDRTLLASILLRIFLHEKLESLLL CTLNDREISMEDEATTLFRATTLASTLMEQYMKATATQFVHHALKDSILKIMESKQSCELSPSKLEKNED VNTNLTHLLNILSELVEKIFMASEILPPTLRYIYGCLQKSVQHKWPTNTTMRTRVVSGFVFLRLICPAIL NPRMFNIISDSPSPIAARTLILVAKSVQNLANLVEFGAKEPYMEGVNPFIKSNKHRMIMFLDELGNVPEL PDTTEHSRTDLSRDLAALHEICVAHSDELRTLSNERGAQQHVLKKLLAITELLQQKQNQYTKTNDVR

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDK**V**

Tag: C-Myc/DDK
Predicted MW: 116.2 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol





RASA1 (NM_002890) Human Recombinant Protein - TP305473M

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: <u>NP 002881</u>

Locus ID: 5921

UniProt ID: <u>P20936</u>, <u>Q59GK3</u>

RefSeq Size: 4402 Cytogenetics: 5q14.3 RefSeq ORF: 3141

Synonyms: CM-AVM; CMAVM; CMAVM1; GAP; p120; p120GAP; p120RASGAP; PKWS; RASA; RASGAP

Summary: The protein encoded by this gene is located in the cytoplasm and is part of the GAP1 family of

GTPase-activating proteins. The gene product stimulates the GTPase activity of normal RAS p21 but not its oncogenic counterpart. Acting as a suppressor of RAS function, the protein enhances the weak intrinsic GTPase activity of RAS proteins resulting in the inactive GDP-bound form of RAS, thereby allowing control of cellular proliferation and differentiation. Mutations leading to changes in the binding sites of either protein are associated with basal cell carcinomas. Mutations also have been associated with hereditary capillary malformations (CM) with or without arteriovenous malformations (AVM) and Parkes Weber syndrome. Alternative splicing results in two isoforms where the shorter isoform, lacking the N-terminal hydrophobic region but retaining the same activity, appears to be abundantly expressed in


placental but not adult tissues. [provided by RefSeq, May 2012]

Protein Families: Druggable Genome

Protein Pathways: Axon guidance, MAPK signaling pathway



Product images:



Coomassie blue staining of purified RASA1 protein (Cat# [TP305473]). The protein was produced from HEK293T cells transfected with RASA1 cDNA clone (Cat# [RC205473]) using MegaTran 2.0 (Cat# [TT210002]).