

Product datasheet for **TP305495L**

WBP2NL (NM_152613) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human WBP2 N-terminal like (WBP2NL), 1 mg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC205495 protein sequence
Red=Cloning site **Green**=Tags(s)

MAVNQSHTENRRGALIPNGESLLKRSPNVELSFPQRSEGSNVFSGRKTGTLFLTSYRVIFITSCSISDPM
LSFMMPFDLMTNLTVEQPVFAANFIKGTIQAAPYGGWEGQATFKLVFRNGGAIEFAQLMVKAASAAARGF
PLRTLNDWFSSMGIYVITGEGNMCTPQMPCSVIVYGAPPAGYGAPPPGYGAPPAGYGAQPVGNEGPPVGY
RASPVRYGAPPLGYGAPPAGYGAPPLGYGAPPLGYGTPLGYGAPPLGYGAPPAGNEGPPAGYRASPAGS
GARPHESTAAQAPENEASLPSASSSQVHS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 31.7 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_689826](#)


Locus ID: 164684



[View online »](#)

UniProt ID:	Q6ICG8
RefSeq Size:	2330
Cytogenetics:	22q13.2
RefSeq ORF:	927
Synonyms:	GRAMD7; PAWP
Summary:	WBP2NL is a sperm-specific WW domain-binding protein that promotes meiotic resumption and pronuclear development during oocyte fertilization (Wu et al., 2007 [PubMed 17289678]). [supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified WBP2NL protein (Cat# [TP305495]). The protein was produced from HEK293T cells transfected with WBP2NL cDNA clone (Cat# [RC205495]) using MegaTran 2.0 (Cat# [TT210002]).