

Product datasheet for **TP305573**

ZSWIM3 (NM_080752) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger, SWIM-type containing 3 (ZSWIM3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205573 protein sequence Red =Cloning site Green =Tags(s)
	MELGSCFKTYEDFKECF SAYKRENRC SFILRDCVSVRFHNLNHGTSIREDILYVQVKFVCIRTQSNRKRT READMCPAYLLLRYNERLDRLFISELNTQHIHGDSKVASPGD TTGKSQKTMCLQRLQPVQPTTKKDLDT AEKSLVEPSFCLDKVQVSSKPEQEGITPSDLAKIAKVMKNFLKVDEGSMASFSVGD SQHLDRLSFQSSKM TDLFIRFPENLLLHRVENTQGHILYAFLVENKERESRVVHFAVLKAETVTSVAKMLSIFTEFNSDWPVKV VWFVDPSFH YRAILQEIFPAARILLSIYHTTRLLEKKLHRSSANPSFKRLMKEALREAVFVTSEASLKNL CQMSQAVLDEDLFNFLQAHWFTCELLWYMHVRKGLLACNTYMDSLDIVTSKVSSLFREQQSLLDCILCFV DYIDFFNTKGLKNLPTPPP KLRARPAS MPLKSKKAFGICGESLTS LPAEETK PDAQVQVQQQS QVPPS QVGMLDTLHQSGSELAYKLCHNEWVQNSTHLVDMAGSSVDVQLLED SHQVSKDGCSCSCSFQQWYHLP CRHILALLHTSQQPVG EAMVCRRWQK KYQYLLGPN GELQDRGMVPNTGQPEKQGRNDMIQDLSRELANLL MQTEGPELEERYSTLRKIVDIWAGPSQPSELFQQPGDFKDVGRLPFLWGKQEEGEGFPATAVMHY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	79.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_542790</u>
Locus ID:	140831
UniProt ID:	<u>Q96MP5</u>
RefSeq Size:	2831
Cytogenetics:	20q13.12
RefSeq ORF:	2088
Synonyms:	C20orf164; PPP1R174

Product images:

Coomassie blue staining of purified ZSWIM3 protein (Cat# TP305573). The protein was produced from HEK293T cells transfected with ZSWIM3 cDNA clone (Cat# [RC205573]) using MegaTran 2.0 (Cat# [TT210002]).