

Product datasheet for **TP305574M**

RBM46 (NM_144979) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RNA binding motif protein 46 (RBM46), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205574 protein sequence Red =Cloning site Green =Tags(s)

MNEENIDGTNGCSKVRTGIQNEAALLALMEKTGYNMVQENGQRKFGGPPPGWEGPPPPRGCEVFGKIPR
DMYEDELVPVFERAGKIYEFRLMMEFSGENRGYAFVMTTKEEAQLAIRILNNYEIRPGKFIGVCSLDN
CRLFIGAIPKEKKKEEILDEMKKVTEGVVDVIVYPSATDKTKNRGFAFVEYESHRAAAMARRKLIPGTGQ
LWGHTIQVDWADPEKEVDEETMQRVKVLYVRNLMISTTEETIKAEFNFKPGAVERVKLRDYAFVHFFN
REDAVAAMSVMNGKCIDGASIEVTLAKPVNKENTWRQHLNGQISPENSENLIVFANKEESHPKTLGKLP
PARLNGQHSPSPPEVERCTYPFYPGTKLTPISMYSLSKNHFNSAVMHLDYCNKNNWAPPEYYLSTTSQ
DGKVLVYKIVIPAIANGSQSYFMPDKLCTTLEDAKELAAQFTLLHLDYFNHRSSINSLSPVSATLSSGT
PSVLPYTSRPYSYPGYPLSPTISLANGSHVGRCLISNQASFF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	59.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_659416](#)

Locus ID: 166863

UniProt ID: [Q8TBY0](#)


RefSeq Size: 2539

Cytogenetics: 4q32.1

RefSeq ORF: 1599

Synonyms: CT68

Product images:



Coomassie blue staining of purified RBM46 protein (Cat# [TP305574]). The protein was produced from HEK293T cells transfected with RBM46 cDNA clone (Cat# [RC205574]) using MegaTran 2.0 (Cat# [TT210002]).