

Product datasheet for **TP305619**

Annexin V (ANXA5) (NM_001154) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human annexin A5 (ANXA5), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205619 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAQVLRGTVTDFPGFDERADAETLRKAMKGLGTDEESILTLTSRSNAQRQEISAAFKTLFGRDLLDDLK SELTGKFEKLIVALKMPSRLYDAYELKHALKGAGTNEKVLTEIIASRTPEELRAIKQVYEEYGSSEDD VVGDTSGYYQRMLVLLQANRDPDAGIDEAQVEQDAQALFQAGELKWGTDEEKFITIFGTRSVSHLRKVF DKYMTISGFQIEETIDRETSGNLEQLLLAVVKSIRSIPAYLAETLYYAMKGAGTDDHTLIRVMVSRSETD LFNIRKEFRKNFATSLYSMIKGDTS GDYKKALLLLCGEDD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	35.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001145</u>
Locus ID:	308



[View online »](#)

UniProt ID: [P08758](#)

RefSeq Size: 1624


Cytogenetics: 4q27

RefSeq ORF: 960

Synonyms: ANX5; ENX2; HEL-S-7; PP4; RPRGL3

Summary: The Annexin 5 gene spans 29 kb containing 13 exons, and encodes a single transcript of approximately 1.6 kb and a protein product with a molecular weight of about 35 kDa. The protein encoded by this gene belongs to the annexin family of calcium-dependent phospholipid binding proteins some of which have been implicated in membrane-related events along exocytotic and endocytotic pathways. Annexin 5 is a phospholipase A2 and protein kinase C inhibitory protein with calcium channel activity and a potential role in cellular signal transduction, inflammation, growth and differentiation. Annexin 5 has also been described as placental anticoagulant protein I, vascular anticoagulant-alpha, endonexin II, lipocortin V, placental protein 4 and anchorin CII. Polymorphisms in this gene have been implicated in various obstetric complications. [provided by RefSeq, Dec 2019]

Product images:



Coomassie blue staining of purified ANXA5 protein (Cat# TP305619). The protein was produced from HEK293T cells transfected with ANXA5 cDNA clone (Cat# [RC205619]) using MegaTran 2.0 (Cat# [TT210002]).