

Product datasheet for TP305660

CYRIB (NM_016623) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 49, member B (FAM49B), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205660 protein sequence Red =Cloning site Green =Tags(s)

MGNLLKVLCTDLEQGPFFLDFENAQPTSESEKEIYNQVNVVLKDAEGILEDLQSYRGAGHEIREAIQHP
ADEKLQEKA WGAVVPLVGLKFKFYEQSRLEAALRGLLGALTSTPYSPTQHLEREQALAKQFAEILHFTL
RFDELKMTNPAIQNDFSYRRTL SRMRINNVP AEGENEVNNELANRMSLFYAEATPMLKTLSDATTKFVS
ENKNLPIENTTDCLSTMASVCRVMLETPEYRSRFTNEETVSFCLRVMVGVIIYDHSVHPVGAFAKTSKID
MKGCIKVLKDQPPNSVEGLLNALRYTTKHLNDETTSKQIKSMLQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	36.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_057707</u>
Locus ID:	51571



[View online »](#)

UniProt ID: [Q9NUQ9](#), [A0A024R9G4](#)

RefSeq Size: 3838


Cytogenetics: 8q24.21

RefSeq ORF: 972

Synonyms: BM-009; CYRI; CYRI-B; FAM49B; L1

Summary: Negatively regulates RAC1 signaling and RAC1-driven cytoskeletal remodeling (PubMed:31285585, PubMed:30250061). Regulates chemotaxis, cell migration and epithelial polarization by controlling the polarity, plasticity, duration and extent of protrusions. Limits Rac1 mediated activation of the Scar/WAVE complex, focuses protrusion signals and regulates pseudopod complexity by inhibiting Scar/WAVE-induced actin polymerization (PubMed:30250061). Protects against Salmonella bacterial infection. Attenuates processes such as macropinocytosis, phagocytosis and cell migration and restrict sopE-mediated bacterial entry (PubMed:31285585). Restricts also infection mediated by Mycobacterium tuberculosis and Listeria monocytogenes (By similarity). Involved in the regulation of mitochondrial dynamics and oxidative stress (PubMed:29059164).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified FAM49B protein (Cat# TP305660). The protein was produced from HEK293T cells transfected with FAM49B cDNA clone (Cat# [RC205660]) using MegaTran 2.0 (Cat# [TT210002]).