

Product datasheet for TP305708

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Adipose Triglyceride Lipase (PNPLA2) (NM_020376) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human patatin-like phospholipase domain containing 2 (PNPLA2), 20

μg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC205708 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MFPREKTWNISFAGCGFLGVYYVGVASCLREHAPFLVANATHIYGASAGALTATALVTGVCLGEAGAKFI EVSKEARKRFLGPLHPSFNLVKIIRSFLLKVLPADSHEHASGRLGISLTRVSDGENVIISHFNSKDELIQ ANVCSGFIPVYCGLIPPSLQGVRYVDGGISDNLPLYELKNTITVSPFSGESDICPQDSSTNIHELRVTNT SIQFNLRNLYRLSKALFPPEPLVLREMCKQGYRDGLRFLQRNGLLNRPNPLLALPPARPHGPEDKDQAVE SAQAEDYSQLPGEDHILEHLPARLNEALLEACVEPTDLLTTLSNMLPVRLATAMMVPYTLPLESALSFTI RLLEWLPDVPEDIRWMKEQTGSICQYLVMRAKRKLGRHLPSRLPEQVELRRVQSLPSVPLSCAAYREALP GWMRNNLSLGDALAKWEECQRQLLLGLFCTNVAFPPEALRMRAPADPAPADPASPQHQLAGPAPLLST

PAPEARPVIGALGL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 55.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.



Adipose Triglyceride Lipase (PNPLA2) (NM_020376) Human Recombinant Protein - TP305708

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 065109

57104 Locus ID: **UniProt ID:** Q96AD5 RefSeq Size: 2443

Cytogenetics: 11p15.5 RefSeq ORF: 1512

Synonyms: 1110001C14Rik; ATGL; FP17548; iPLA2zeta; PEDF-R; TTS-2.2; TTS2


Summary: This gene encodes an enzyme which catalyzes the first step in the hydrolysis of triglycerides in

adipose tissue. Mutations in this gene are associated with neutral lipid storage disease with

myopathy. [provided by RefSeq, Jul 2010]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified PNPLA2 protein (Cat# TP305708). The protein was produced from HEK293T cells transfected with PNPLA2 cDNA clone (Cat# [RC205708]) using MegaTran 2.0 (Cat# [TT210002]).