

## Product datasheet for TP305732

### C2orf15 (NM\_144706) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 2 open reading frame 15 (C2orf15), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205732 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MSSTLGKLSNQEETLPLLKKVPANYFHICSAILMGFSLSKSATQVSAIHMSKVDHDLIRGTEKSRLEP ATQLFQNTKKIRLEDTNQENFTRIEGTGTGSLSGKALGSVYVKESDGLEMTDVE
	<b>TR</b> TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_653307</a>
Locus ID:	150590
UniProt ID:	<a href="#">Q8WU43</a> , <a href="#">B2RAC6</a>
RefSeq Size:	1430
Cytogenetics:	2q11.2



[View online »](#)

RefSeq ORF: 375

### Product images:



Coomassie blue staining of purified C2orf15 protein (Cat# TP305732). The protein was produced from HEK293T cells transfected with C2orf15 cDNA clone (Cat# [RC205732]) using MegaTran 2.0 (Cat# [TT210002]).