

Product datasheet for **TP305764**

FBLIM1 (NM_017556) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human filamin binding LIM protein 1 (FBLIM1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205764 protein sequence Red =Cloning site Green =Tags(s)

MASKPEKRVASSVFITLAPPRRDVAVAEEVRQAVCEARRGRPWEAPAPMKTPEAGLAGRPSWTTTPGRAA
ATVPAAPMQLFNGGCPPPPVLDGEDVLPDLLPPPPPPVLLPSEEEAPAPMGASLIADLEQLHLSP
PPPPQAPAEGPSVQPGPLRPMEELPPPPAEPVEKGASTDICAFCFKTVFPRELAVEAMKRQYHAQCFT
CRTCRRQLAGQSFYQKDGRLCEPCYQDTLERCCKGKCGEVWRDHIIRALGQAFHPSCFTCVTCARCIGDES
FALGSQNEVYCLDDFYRKFAPVCSICENPIIPRDGKDAFKIECMGRNFHENCYRCEDCRILLSVEPTDQG
CYPLNNHLFCKPCHVKRSAAGCC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	40.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_060026</u>



[View online »](#)

Locus ID:	54751
UniProt ID:	Q8WUP2
RefSeq Size:	3363
Cytogenetics:	1p36.21
RefSeq ORF:	1119
Synonyms:	CAL; FBLP-1; FBLP1

Summary: This gene encodes a protein with an N-terminal filamin-binding domain, a central proline-rich domain, and, multiple C-terminal LIM domains. This protein localizes at cell junctions and may link cell adhesion structures to the actin cytoskeleton. This protein may be involved in the assembly and stabilization of actin-filaments and likely plays a role in modulating cell adhesion, cell morphology and cell motility. This protein also localizes to the nucleus and may affect cardiomyocyte differentiation after binding with the CSX/NKX2-5 transcription factor. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified FBLIM1 protein (Cat# TP305764). The protein was produced from HEK293T cells transfected with FBLIM1 cDNA clone (Cat# [RC205764]) using MegaTran 2.0 (Cat# [TT210002]).