

## Product datasheet for **TP305771L**

### PGM1 (NM\_002633) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphoglucomutase 1 (PGM1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205771 representing NM_002633 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MVKIVTVKTQAYQDQKPGTSGLRKRKVKVFQSSANYAENFIQSIISTVEPAQRQEATLVVGGDGRFYMKEA IQLIARIAAANGIGRLVIGQNGILSTPAVSCIIRKKAIGGIILTASHNPGGPNDFGIKFNISNGGPAP EAITDKIFQISKTIIEYAVCPDLKVDLGLVGLKQKQFDLENKFKPFTVEIVDSVEAYATMLRSIFDFSALKE LLSGPNRLKIRIDAMHGVPYVKKILCEELGAPANSVNCVPLEDFGGHHPDPNLTYAADLVETMKSGE HDFGAAFDGDGDRNMILGKHGFFVNPDSVAVIAANIFSIPYFQQTGVRGFFARSMPTSGALDRVASATKI ALYETPTGWKFFGNLMDASKLSLCGEESFGTGSDHIREKDGLWAVLAWLSILATRKQSVEDILKDHQKY GRNFFTRYDYEEVEAEGANKMMKDLEALMFDRSFVGKQFSANDKVYTVKADNFEYSDPVDGSIQRNQL RLIFTDGSRIVFRLSGTGSAGATIRLYIDSYEKDVAKINQDPQVMLAPLISIALKVSQQLQERTGRTAPT IT</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	61.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_002624</a>
<b>Locus ID:</b>	5236
<b>UniProt ID:</b>	<a href="#">P36871</a>
<b>RefSeq Size:</b>	2487
<b>Cytogenetics:</b>	1p31.3
<b>RefSeq ORF:</b>	1686
<b>Synonyms:</b>	CDG1T; GSD14
<b>Summary:</b>	The protein encoded by this gene is an isozyme of phosphoglucomutase (PGM) and belongs to the phosphohexose mutase family. There are several PGM isozymes, which are encoded by different genes and catalyze the transfer of phosphate between the 1 and 6 positions of glucose. In most cell types, this PGM isozyme is predominant, representing about 90% of total PGM activity. In red cells, PGM2 is a major isozyme. This gene is highly polymorphic. Mutations in this gene cause glycogen storage disease type 14. Alternativley spliced transcript variants encoding different isoforms have been identified in this gene.[provided by RefSeq, Mar 2010]
<b>Protein Pathways:</b>	Amino sugar and nucleotide sugar metabolism, Galactose metabolism, Glycolysis / Gluconeogenesis, Metabolic pathways, Pentose phosphate pathway, Starch and sucrose metabolism

### Product images:



Coomassie blue staining of purified PGM1 protein (Cat# [TP305771]). The protein was produced from HEK293T cells transfected with PGM1 cDNA clone (Cat# [RC205771]) using MegaTran 2.0 (Cat# [TT210002]).