

## Product datasheet for **TP305792L**

### ABHD12B (NM\_181533) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human abhydrolase domain containing 12B (ABHD12B), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205792 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MLGIWHTVPSCRGEDAKGKDCCWYEALRDGNPIIVYLHGSAEHRAASHRLKLVKVLSDGGFHVLSVDYR  
GFGDSTGKPTTEGLTTDAICVYEWTKARSGITPVCLWGHSLGTGVATNAAKVLEEKGCPVDAIVLEAPFT  
NMWVASINYPLLKIYRNIPGFLRTLMDALRKDKIIFPNDENVKFLSSPLLILHGEDDRTVPLEYGKKLYE  
IARNAYRNKERVKMVFPPGFQHNLLCKSPTLLITVRDFLSKQWS

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


Tag:	C-Myc/DDK
Predicted MW:	28.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_853511</a></u>
Locus ID:	145447



[View online »](#)

UniProt ID: [Q7Z5M8](#)  
RefSeq Size: 1580  
Cytogenetics: 14q22.1  
RefSeq ORF: 765  
Synonyms: BEM46L3; C14orf29; c14\_5314

### Product images:



Coomassie blue staining of purified ABHD12B protein (Cat# [TP305792]). The protein was produced from HEK293T cells transfected with ABHD12B cDNA clone (Cat# [RC205792]) using MegaTran 2.0 (Cat# [TT210002]).