

Product datasheet for **TP305833M**

CYB5R1 (NM_016243) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cytochrome b5 reductase 1 (CYB5R1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205833 protein sequence Red =Cloning site Green =Tags(s)
	 MGIQTSPVLLASLGVGLVTLGLAVGSYLVRSSRRPQVTLDPNEKYLLRLLDKTTVSHNTRKFRFALPT AHHTLGLPVGKHIYLSTRIDGSLVIRPYTPVTSDEDQGYVDLVIKVYLKGVHPKPFEGGKMSQYLDLTKV GDVVEFRGPSGLLTYTGKGFNIQPNKKSPPEPRVAKKLGMIAGGTGITPMLQLIRAILKVPEDPTQCFL LFANQTEKDIILREDLEELQARYPNRFLWFTLDHPPKDWAYSKGFVTADMIREHLPAGDDVLVLLCGP PPMVQLACHPNLCLKGYSQKMRFTY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	33.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_057327</u>
Locus ID:	51706



[View online »](#)

UniProt ID:	<u>Q9UHQ9</u>
RefSeq Size:	1675
Cytogenetics:	1q32.1
RefSeq ORF:	915
Synonyms:	B5R.1; B5R1; B5R2; humb5R2; NQO3A2
Summary:	NADH-cytochrome b5 reductases are involved in desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction. [UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Amino sugar and nucleotide sugar metabolism

Product images:



Coomassie blue staining of purified CYB5R1 protein (Cat# [TP305833]). The protein was produced from HEK293T cells transfected with CYB5R1 cDNA clone (Cat# [RC205833]) using MegaTran 2.0 (Cat# [TT210002]).