

Product datasheet for **TP305847**

SNX8 (NM_013321) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human sorting nexin 8 (SNX8), 20 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC205847 protein sequence
Red=Cloning site **Green**=Tags(s)

MTGRAMDPLPAAAVGAAAEAEDEEADPPASDLPTPQAIEPQAIQVQVPAPSRMQMPQGNPLLSHTLQE
LLARDTVQVELIPEKKGLFLKHVEYEVSSQRFKSSVYRRYNDFFVVFQEMLLHKFPYRMVVPALPPKRMGLGA
DREFIEARRRALKRFVNLVARHPLFSEDVVLKLFSLFSGSDVQNKLKESAQCVGDEFNCKLATRAKDFL
PADIQAQFAISRELIRNIYNSFHKLDRDRAERIASRAIDNAAADLLIFGKELSAIGSDTTPPLPSWAALNSST
WGLKQALKGLSVEFALLADKAAQQKGQEENDVVEKLNFLDLLQSYKDLKERHEKGVLHKHQRALHKYS
LMKRQMMSATAQNREPESVEQLESRIVEQENAIQTMELRNYFSLYCLHQETQLIHVYLPLTSHILRAFVN
SQIQGHKEMSKVWNDLRPKLSCLFAGPHSTLTTPCSPPEDGLCPH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 52.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


RefSeq: [NP_037453](#)



[View online »](#)

Locus ID:	29886
UniProt ID:	Q9Y5X2
RefSeq Size:	4756
Cytogenetics:	7p22.3
RefSeq ORF:	1395
Synonyms:	Mvp1
Summary:	May be involved in several stages of intracellular trafficking. May play a role in intracellular protein transport from early endosomes to the trans-Golgi network.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified SNX8 protein (Cat# TP305847). The protein was produced from HEK293T cells transfected with SNX8 cDNA clone (Cat# [RC205847]) using MegaTran 2.0 (Cat# [TT210002]).