

## Product datasheet for **TP305880M**

### KDEL Receptor (KDEL1) (NM\_006801) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human KDEL (Lys-Asp-Glu-Leu) endoplasmic reticulum protein retention receptor 1 (KDEL1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205880 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MNLFRFLGDLSHLLAIILLLLKIWKSRSCAGISGKSQVLFAVFTARYLDLFTNYISLYNTCMKWYIAC SFTTVWLIYSKFKATYDGNHDTFRVEFLVPTAILAFLVNHDFPLEILWTFESIYLESVAILPQLFMVSK TGEAETITSHYLFALGVYRTLFLFNWIWRYHFEGFFDLIAIVAGLVQTVLYCDDFFLYITKVLKGGKLSL PA  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	24.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_006792</a></u>
Locus ID:	10945



[View online »](#)

UniProt ID: [P24390](#)

RefSeq Size: 1575

Cytogenetics: 19q13.33

RefSeq ORF: 636


Synonyms: ERD2; ERD2.1; HDEL; PM23

**Summary:** Retention of resident soluble proteins in the lumen of the endoplasmic reticulum (ER) is achieved in both yeast and animal cells by their continual retrieval from the cis-Golgi, or a pre-Golgi compartment. Sorting of these proteins is dependent on a C-terminal tetrapeptide signal, usually lys-asp-glu-leu (KDEL) in animal cells, and his-asp-glu-leu (HDEL) in *S. cerevisiae*. This process is mediated by a receptor that recognizes, and binds the tetrapeptide-containing protein, and returns it to the ER. In yeast, the sorting receptor encoded by a single gene, ERD2, which is a seven-transmembrane protein. Unlike yeast, several human homologs of the ERD2 gene, constituting the KDEL receptor gene family, have been described. The protein encoded by this gene was the first member of the family to be identified, and it encodes a protein structurally and functionally similar to the yeast ERD2 gene product. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Vibrio cholerae infection

### Product images:



Coomassie blue staining of purified KDEL1 protein (Cat# [TP305880]). The protein was produced from HEK293T cells transfected with KDEL1 cDNA clone (Cat# [RC205880]) using MegaTran 2.0 (Cat# [TT210002]).