

Product datasheet for **TP305890**

Aromatase (CYP19A1) (NM_000103) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human cytochrome P450, family 19, subfamily A, polypeptide 1 (CYP19A1), transcript variant 1, full length, with C-terminal MYC/DDK tag, expressed in HEK293 cells, 20 µg
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	>RC205890 protein sequence Red =Cloning site Green =Tags(s)
	<p>MVLEMLNPIHYNITSIVPEAMPAATMPVLLLTLGLFLLVWNYEGTSSIPGPGYCMGIGPLISHGRFLWMGI GSACNYNRVYGEFMRVWISGEETLIISKSSSMFHIMKHNHYSSRFGSKLGLQCIGMHEKGIIFNNNP WKTTRPFFMKALSGPGLVRMVTVCAESLKTHLDRLEEVTNESGYVDVLTLLRRVMLDTSNTLFLRISLDE SAIVVKIQGYFDAWQALLIKPDIFFKISWLYKKEYKSVKDLKDAIEVLIAEKRRRISTEEKLEECMDFAT ELILAEKRGDLTRENVNQCILEMLIAAPDTMSVSLFFMLFLIAKHPNVEEAIKEIQTVIGERDIKIDDI QKLKVMENFIYESMRYQPVVDLVMRKALEDDVIDGYPVKKGTNIIINIGRMHRLEFFPKPNEFTLENFAK NVPYRYFQPFQFGPRGCAGKYIAMVMMKAILVTLRRFHVKTLLQGCVESIQKIHDLSLHPDETKNMLEM IFTPRNSDRCLEH</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	57.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_000094</u>
Locus ID:	1588
UniProt ID:	<u>P11511, A0A024R5S8, Q8TCA4, Q8IYG4</u>
RefSeq Size:	4422
Cytogenetics:	15q21.2
RefSeq ORF:	1509
Synonyms:	ARO; ARO1; CPV1; CYAR; CYP19; CYPXIX; P-450AROM
Summary:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and catalyzes the last steps of estrogen biosynthesis. Mutations in this gene can result in either increased or decreased aromatase activity; the associated phenotypes suggest that estrogen functions both as a sex steroid hormone and in growth or differentiation. Alternative promoter use and alternative splicing results in multiple transcript variants that have different tissue specificities. [provided by RefSeq, Dec 2016]
Protein Families:	Druggable Genome, P450
Protein Pathways:	Androgen and estrogen metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified CYP19A1 protein (Cat# TP305890). The protein was produced from HEK293T cells transfected with CYP19A1 cDNA clone (Cat# [RC205890]) using MegaTran 2.0 (Cat# [TT210002]).