

Product datasheet for TP305976

PIP4K2C (NM_024779) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphatidylinositol-5-phosphate 4-kinase, type II, gamma (PIP4K2C), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205976 protein sequence Red=Cloning site Green=Tags(s)

MASSSVPPATVSAATAGPGPGFGFASKTKKKHFVQQKVKVFRAADPLVGVFLWGVVAHSINELSQVPPPVM
LLPDDFKASSKIKVNNHLFHRENLP SHFKFKEYCPQVFRNLRDRFGIDDQDYLVSLTRNPPSESEGSDGR
FLISYDRTLVIKEVSSEDIADMHSNLSNYHQYIVKCHGNTLLPQFLGMYRVSDNEDSYMLVMRNMFSHR
LPVHRKYDLKGSLSREASDKEKVKELPTLRDMDFLNKNQKVYIGEEKKIFLEKLKRDVEFLVQLKIMD
YSLLLGIHDIIRGSEPEEEAPVREDESEVDGDCSLTGPPALVGSYGTSPEGIGGGYIHSRPLGPGEFESF
IDVYAIRSAEGAPQKEVYFMGLIDILTQYDAKKKAAHAAKTVKHGAGAEISTVHPEQYAKRFLDFITNIF
A

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	47.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	NP_079055
Locus ID:	79837
UniProt ID:	Q8TBX8
RefSeq Size:	3229
Cytogenetics:	12q13.3
RefSeq ORF:	1263
Synonyms:	PIP5K2C
Summary:	May play an important role in the production of Phosphatidylinositol bisphosphate (PIP2), in the endoplasmic reticulum.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome
Protein Pathways:	Inositol phosphate metabolism, Phosphatidylinositol signaling system, Regulation of actin cytoskeleton

Product images:



Coomassie blue staining of purified PIP4K2C protein (Cat# TP305976). The protein was produced from HEK293T cells transfected with PIP4K2C cDNA clone (Cat# [RC205976]) using MegaTran 2.0 (Cat# [TT210002]).