

Product datasheet for **TP305986M**

Gephyrin (GPHN) (NM_020806) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human gephyrin (GPHN), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC205986 protein sequence Red =Cloning site Green =Tags(s)

MATEGMILTNDHQIRVGVLTVSDSCFRNLAEDRSGINLKDLVQDPSLLGGTISAYKIVPDEIEEIKETL
IDWCDEKELNLILTTGGTGAFPRDVTPEATKEVIEREAPGMALAMLMGSLNVTPLGMLSRPVCGIRGKTL
IINLPGSKKGSQECFQFILPALPHAILLRDAIVKVKVHDELEDLPSPPPLSPPTTSPHKQTEDKGV
QCEEEEEKKGSGVASTEDSSSSHITAAIAAKKHFPYTSAPVMAHGEQPIPLINYSHHSTDERIPDS
IISRGVQVLPRTASLSTTPSESPRAQATSRLSTASCPTPKVQSRCSSKENILRASHSAVDITKVARHR
MSPFPLTSMDFKAFITVLEMTPLVLTGTEIINYRDGMGRVLAQDVYAKDNLPPFPASVKDGYAVRAADGPGDR
FIIGESQAGEQPTQTVMPGQVMRVTTGAPIPCGADAVVQVEDTELIRESDDGTEELEVRILVQARPGQDI
RPIGHDIKRGECVLAAGTHMGPSEIGLLATVGVTEVEVNFVWAVMSTGNELNPPEDDLLPGKIRDSNR
STLLATIQEHGYPTINLGIQVGNPDDLLNALNEGISRADVIITSGGVSMGEKDYKQVLDIDLHAQIHFG
RVFMKPLPTTFATLDIDGVRKIIFALPGNPVSAVTCNLFVWPALRKMVGILDRPTIHKARLSCDVKL
DPRPEYHRCILTWHHQEPLPWAQSTGNQMSSRLMSMRANGLLMLPPKTEQYVELHKGEVVDVMVIGRL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	83.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online >](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_065857
Locus ID:	10243
UniProt ID:	Q9NQX3
RefSeq Size:	4318
Cytogenetics:	14q23.3-q24.1
RefSeq ORF:	2307
Synonyms:	GEPH; GPH; GPHRYN; HKPX1; MOCODC
Summary:	This gene encodes a neuronal assembly protein that anchors inhibitory neurotransmitter receptors to the postsynaptic cytoskeleton via high affinity binding to a receptor subunit domain and tubulin dimers. In nonneuronal tissues, the encoded protein is also required for molybdenum cofactor biosynthesis. Mutations in this gene may be associated with the neurological condition hyperplexia and also lead to molybdenum cofactor deficiency. Numerous alternatively spliced transcript variants encoding different isoforms have been described; however, the full-length nature of all transcript variants is not currently known. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified GPHN protein (Cat# [TP305986]). The protein was produced from HEK293T cells transfected with GPHN cDNA clone (Cat# [RC205986]) using MegaTran 2.0 (Cat# [TT210002]).