

Product datasheet for TP306085

DDX55 (NM_020936) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human DEAD (Asp-Glu-Ala-Asp) box polypeptide 55 (DDX55), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206085 protein sequence Red =Cloning site Green =Tags(s)
	<p>MEHVTEGSWESLPVPLHPQVLGALRELGFYMTVPVQSATIPLFMRNKDVAAEAVTGSGKTLAFVIPILEI LLRREEKLKKSQVGAIITPTRELAIQIDEVLSHFTKHFPEFSQILWIGGRNPGEDVERFKQQGGNIIVA TPGRLEDLFRRKAEGDLASCVRSLDVLVLDEADRLLDMGFEASINTILEFLPKQRRRTGLFSATQTQVE NLVRAGLRNPVRVSVEKEGVAASSAQKTPSRLNYYMVCKADEKFNQLVHFLRNHKQEKHLVFFSTCACV EYYGKTLEVLVKGVKIMCIHGKMKYKRNKIFMEFRKLQSGILVCTDVMARGIDIPEVNWWVLQYDPPSNAS AFVHRCGRTARIGHGGSALVLLPMEESYINFLAINQKCPLQEMKQQRNTADLLPKLKSMALADRAVFEK GMKAFVSYVQAYAKHECNLIFRLKDLDFASLARGFALLRMPKMPQLRGKQFPDFVPVDVNTDTIPFKDKI REKQRQKLLQQRREKTENEGRRKFIKNKAWSKQKAKKEKKKKMNEKRKREEGSDIEDEDMEELLNDTRL LKLLKKGKITEEEFEKGLLTGKRTIKTVDLGISDLEDGC</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	68.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_065987</u>
Locus ID:	57696
UniProt ID:	<u>Q8NHQ9</u>
RefSeq Size:	2638
Cytogenetics:	12q24.31
RefSeq ORF:	1800
Summary:	This gene encodes a member of protein family containing a characteristic Asp-Glu-Ala-Asp (DEAD) motif. These proteins are putative RNA helicases, and may be involved in a range of nuclear processes including translational initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Multiple alternatively spliced transcript variants have been found for this gene. Pseudogenes have been identified on chromosomes 1 and 12. [provided by RefSeq, Feb 2016]

Product images:



Coomassie blue staining of purified DDX55 protein (Cat# TP306085). The protein was produced from HEK293T cells transfected with DDX55 cDNA clone (Cat# [RC206085]) using MegaTran 2.0 (Cat# [TT210002]).