

Product datasheet for **TP306131L**

GPD1L (NM_015141) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glycerol-3-phosphate dehydrogenase 1-like (GPD1L), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206131 protein sequence Red =Cloning site Green =Tags(s)

MAAAPLKVCIVGSGNWGS AVAKIIGNNVKKLQKFASTVKMWWFEETVNGRKLTDIINNDHENVKYLPGHK
LPENWVAMSNLSEAVQDADLLV FVIPHQFIHRICDEITGRVPKKALGITLIKIDEGPEGLKLISDIIRE
KMGIDISVLMGANIANEVAAEKFCETTIGSKVMENGLLFKELLQTPNFRITVDDADTVELCGALKNIVA
VGAGFCDGLRCGDNTKAAVIRLGLMEMIAFARIFCKGQVSTATFLESCGVADLITTCYGGRRNRVAEFA
RTGKTIEELEKEMLNQKLQGPQTSAEVYRILKQKGLLDKFPLFTAVYQICYESRPVQEMLSCLSQSHPEH
T

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	38.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_055956
Locus ID:	23171



[View online »](#)

UniProt ID: [Q8N335](#)

RefSeq Size: 4068

Cytogenetics: 3p22.3


RefSeq ORF: 1053

Synonyms: GPD1-L

Summary: The protein encoded by this gene catalyzes the conversion of sn-glycerol 3-phosphate to glycerone phosphate. The encoded protein is found in the cytoplasm, associated with the plasma membrane, where it binds the sodium channel, voltage-gated, type V, alpha subunit (SCN5A). Defects in this gene are a cause of Brugada syndrome type 2 (BRS2) as well as sudden infant death syndrome (SIDS). [provided by RefSeq, Jul 2010]

Protein Pathways: Glycerophospholipid metabolism

Product images:



Coomassie blue staining of purified GPD1L protein (Cat# [TP306131]). The protein was produced from HEK293T cells transfected with GPD1L cDNA clone (Cat# [RC206131]) using MegaTran 2.0 (Cat# [TT210002]).