

Product datasheet for **TP306155L**

CCDC117 (NM_173510) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human coiled-coil domain containing 117 (CCDC117), 1 mg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC206155 protein sequence
Red=Cloning site **Green**=Tags(s)

MAALGRPFSGPLPLSGGSDFLQPPQPAFPGRAFPFGADGAELAPRPGPRAVPSSPAGSAARGRVSVHCKKK
HKREEEEDDDCPVRKKRITEAELCAGPNDWILCAHQDVEGHGVPNSVSGLSIPGILDVICEEMDQTTGEP
QCEVARRKLQEIEDRIIDEDEEVEADRNVNHLPSLVLSDTMKTGLKREFDEVFTKKMIESMSRPSMELVL
WKPLPELLSDPKPKSSNTKNYTGESQAKHVAAGTAFPQRTELFSEPRPTGMSLYNSLETATSTEEEMEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 30.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_775781](#)


Locus ID: 150275

UniProt ID: [Q8IWD4](#), [A0A024R1C5](#)



[View online »](#)

RefSeq Size: 3991
Cytogenetics: 22q12.1
RefSeq ORF: 837
Synonyms: dj366L4.1

Product images:

Coomassie blue staining of purified CCDC117 protein (Cat# [TP306155]). The protein was produced from HEK293T cells transfected with CCDC117 cDNA clone (Cat# [RC206155]) using MegaTran 2.0 (Cat# [TT210002]).