

Product datasheet for **TP306173**

RAB11FIP5 (NM_015470) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RAB11 family interacting protein 5 (class I) (RAB11FIP5), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206173 representing NM_015470 Red =Cloning site Green =Tags(s)

MALVRGAEPAAGPSRWLPTHVQVTVLRARGLRGKSSGAGSTSDAYTVIQVGREKYSTSVVEKTHGCPEWR
EECSFELPPGALDGLLRAQEADAGPAPWAASSAAACELVLTMMHRSLLIGVDKFLGQATVALDEVFGAGRA
QHTQWYKHLHSPGKKEKERGEIEVTIQFTRNNLSAMFDLSMKDKPRSPFSKIRDKMKGKKKYDLESASAI
LPSSAIEDPDLGSLGKMGKAKGFFLRNKLKSSLTQSNTSLGSDSTLSSASGLAYQGPGAELLTRSPSR
SSWLSTEGGRDSAQSPKLFTHKRTYSDEANQMRVAPPRALLDLQGHLDAAASRSLCVNGSHIYNEEPQGP
VRHRSSISGSLPSSGSLQAVSSRFSEEGPRSTDDTWPRGSRNSNSSEAVLGQEELSAQAKVLAPGASHPG
EEEGARLPEGKPVQVATPIVASSEVAEKEGARKEERKPRMGLFHHHHQGLSRSELGRRSSLGEKGGPIL
GASPHHSSGEGKAKSSWFLREAKDPTQKPSHPVKPLSAAPVEGSPDRKQSRSSLSIALSSGLEKLT
VTSGSIQPVTQAPQAGQMVDTKRLKDSAVLDQSAKYNNLTHDELISLLLQRERELSQRDEHVQELESYID
RLLVRIMETSPTLLQIPPGPPK


SGPTRRRLE**QKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	70.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online >](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_056285
Locus ID:	26056
UniProt ID:	Q9BXF6 , Q9UFM0
RefSeq Size:	4340
Cytogenetics:	2p13.2
RefSeq ORF:	1956
Synonyms:	gaf-1; GAF1; pp75; rab11-FIP5; RIP11
Summary:	Rab effector involved in protein trafficking from apical recycling endosomes to the apical plasma membrane. Involved in insulin granule exocytosis. May regulate V-ATPase intracellular transport in response to extracellular acidosis.[UniProtKB/Swiss-Prot Function]
Protein Pathways:	Endocytosis

Product images:

Coomassie blue staining of purified RAB11FIP5 protein (Cat# TP306173). The protein was produced from HEK293T cells transfected with RAB11FIP5 cDNA clone (Cat# [RC206173]) using MegaTran 2.0 (Cat# [TT210002]).