

Product datasheet for **TP306195M**

ANKRD29 (NM_173505) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ankyrin repeat domain 29 (ANKRD29), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206195 protein sequence Red =Cloning site Green =Tags(s)
	MCRMSFKKETPLANAFAWAARRGNLALLKLLNSGRVDVDCRDSHGTTLLMVAAYAGHIDCVRELVLQG A DINLQRESGTTALFFAAQQGHNDVVRFLFGFGASTEFR TKDEGTALLAASQYGHMQVETLLKHGANIHD QLYDGATALFLAAQGGYLDVIRLLLASGAKVNQPRQDGTAPLWIASQMGHSEVVRVMLLRGADRDAARN D GTTALLKAANKGYNDVIKELLKFSPTLGILKNGTSALHAAVLSGNIKTVALLEAGADPSLRNKANELPA ELTKNERILRLRSKEGPRKS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	32.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_775776</u>



[View online »](#)

Locus ID: 147463
 UniProt ID: [Q8N6D5](#)
 RefSeq Size: 2291
 Cytogenetics: 18q11.2
 RefSeq ORF: 903

Product images:


Coomassie blue staining of purified ANKRD29 protein (Cat# [TP306195]). The protein was produced from HEK293T cells transfected with ANKRD29 cDNA clone (Cat# [RC206195]) using MegaTran 2.0 (Cat# [TT210002]).