

Product datasheet for TP306197L

LAYN (NM_178834) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human layilin (LAYN), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206197 protein sequence Red=Cloning site Green=Tags(s)
	MRPGTALQAVLLAVLLVGLRAATGRLLSGQPVCRRGGTQRPCYKVIYFHDTSRRLNFEEAKEACRRDGGQL VSIESEDEQKLIKFIENLLPSDGFWIGLRRREEKQSNSTACQDLYAWTDGSISQFRNWVYVDEPSCGSE VCVVMYHQPSAPAGIGGPYMFQWDDRCNMKNNFICKYSDEKPAVPSREAEGETELTTPVLPEETQEED AKKTFKESREAALNLAYILIPSIPLLLLLLVTTVVCWWVICRKRKREQDPSTKKQHTIWSPHQGNP LEVYNVIRKQSEADLAETRPDLKNISFRVCSGEATPDDMSCDYDNMAVNPSESGFMTLVSVESGFVTNDI YEFSPDQMGRSKESGWVENEIYGY
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	42.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_849156
Locus ID:	143903



[View online »](#)

UniProt ID:	Q6UX15
RefSeq Size:	2748
Cytogenetics:	11q23.1
RefSeq ORF:	1122
Summary:	Receptor for hyaluronate.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified LAYN protein (Cat# [TP306197]). The protein was produced from HEK293T cells transfected with LAYN cDNA clone (Cat# [RC206197]) using MegaTran 2.0 (Cat# [TT210002]).