

Product datasheet for **TP306200**

ZCCHC24 (NM_153367) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger, CCHC domain containing 24 (ZCCHC24), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206200 protein sequence Red =Cloning site Green =Tags(s)
	 MSLLSAIDTSAASVYQPAQLLNWVYLSLQDTHQASAFDAFRPVPTAGAAPPELAFGKGRPEQLGSPLHSS YLNSFFQLQRGEALSNSVYKGASPYGSLNNIADGLSSLTEHFSDLTLTSEARKPSKRPPPNYLCHLCFNK GHYIKDCPQARPKGEGLTPYQGKKRCFGEYKCPKCKRKMMSGNSWANMGQECIKCHINVYPHKQRPLEKP DGLDVSDQSKEHPQHLCEKCKVLGYCRRVQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	26.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_699198
Locus ID:	219654
UniProt ID:	Q8N2G6 , A0A024QZP0



[View online »](#)

RefSeq Size: 4933
Cytogenetics: 10q22.3
RefSeq ORF: 723
Synonyms: C10orf56; Z3CXXC8

Product images:

Coomassie blue staining of purified ZCCHC24 protein (Cat# TP306200). The protein was produced from HEK293T cells transfected with ZCCHC24 cDNA clone (Cat# [RC206200]) using MegaTran 2.0 (Cat# [TT210002]).