

## Product datasheet for **TP306202M**

### **C10orf82 (NM\_144661) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 10 open reading frame 82 (C10orf82), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206202 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MEPSKTFMRNLPITPGYSGFVPLSCQGMSKEDDMNHCVKTFQEKTQRYKEQLRELCCA VATAPKLPVN SEETVLQALHQYNLQYHPLILECKYVKKPLQEPPIPGWAGYLPRAKVTEFGCGTRYVMANKCYKDFLEI TERAKKAHLKPYEE
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	17.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_653262</a>
Locus ID:	143379
UniProt ID:	<a href="#">Q8WW14</a>
RefSeq Size:	1009




[View online »](#)

Cytogenetics: 10q25.3

RefSeq ORF: 462

### Product images:



Coomassie blue staining of purified C10orf82 protein (Cat# [TP306202]). The protein was produced from HEK293T cells transfected with C10orf82 cDNA clone (Cat# [RC206202]) using MegaTran 2.0 (Cat# [TT210002]).