

Product datasheet for **TP306225M**

C5orf35 (SETD9) (NM_153706) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human chromosome 5 open reading frame 35 (C5orf35), 100 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC206225 protein sequence
Red=Cloning site **Green**=Tags(s)

MPGRLLRGLWQRWRRYKYRFVPWIALNLSHNPRTLRYVPEESKDKVISDEDVLGTLTKVFQALFLNDFNK
QSEILTMLPESVKSKYQDLLAVEHQGVKLRENHRHQQSTFKPEEILYKTLGFSVAQATSSLISAGKGVFV
TKGLVPKGAVSMYPGTVYQKYEPIFFQSIGNPFIFRCLDGVLDGNDKGISKVVYRSCNGRDRLGPKM
SDSTWLTSEIHNPPLAVGQYVNNCSNDRAANVCYQEFDVPVAVFPIELKQYLPNIAYSYDKQSPLRCVVLVA
LRDINQGEELFSNYTIVS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 33.9 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_714917](#)
Locus ID: 133383



[View online »](#)

UniProt ID:	<u>Q8NE22</u>
RefSeq Size:	1722
Cytogenetics:	5q11.2
RefSeq ORF:	897
Synonyms:	C5orf35

Product images:

Coomassie blue staining of purified SETD9 protein (Cat# [TP306225]). The protein was produced from HEK293T cells transfected with SETD9 cDNA clone (Cat# [RC206225]) using MegaTran 2.0 (Cat# [TT210002]).