

## Product datasheet for **TP306241L**

### SH3BGRL2 (NM\_031469) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human SH3 domain binding glutamic acid-rich protein like 2 (SH3BGRL2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206241 representing NM_031469 Red=Cloning site Green=Tags(s)
	MVIRVFIASSSGFVAIKKKQDWRFLFLEANKIEFEEVDITMSEEQRQWYKYNVPPEKKPTQGNPLPPQIF NGDRYCGDYDSFFESKESNTVFSFLGLKPRLASKAEP
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	12.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_113657</a>
Locus ID:	83699
UniProt ID:	<a href="#">Q9UJC5</a>
RefSeq Size:	4676




[View online »](#)

Cytogenetics: 6q14.1

RefSeq ORF: 321

### Product images:



Coomassie blue staining of purified SH3BGRL2 protein (Cat# [TP306241]). The protein was produced from HEK293T cells transfected with SH3BGRL2 cDNA clone (Cat# [RC206241]) using MegaTran 2.0 (Cat# [TT210002]).