

## Product datasheet for **TP306249**

### EIF5A2 (NM\_020390) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human eukaryotic translation initiation factor 5A2 (EIF5A2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206249 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MADEIDFTTGDAGASSTYPMQCSALRKNGFVWLKGRPCKIVEMSTSKTGKHGHAKVHLVGDIFTGKKYE DICPSTHNMDVPNIKRNDYQLICIQDGYLSLLTETGEVREDLKLPEGELGKEIEGKYNAGEDVQVSVMCA MSEYAVAIKPKC  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	16.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_065123</a></u>
Locus ID:	56648
UniProt ID:	<u><a href="#">Q9GZV4</a></u>
RefSeq Size:	5537



[View online »](#)


Cytogenetics: 3q26.2

RefSeq ORF: 459

Synonyms: EIF-5A2; eIF5AII

**Summary:** mRNA-binding protein involved in translation elongation. Has an important function at the level of mRNA turnover, probably acting downstream of decapping. Involved in actin dynamics and cell cycle progression, mRNA decay and probably in a pathway involved in stress response and maintenance of cell wall integrity. Functions as a regulator of apoptosis. Mediates effects of polyamines on neuronal process extension and survival. May play an important role in brain development and function, and in skeletal muscle stem cell differentiation (By similarity).[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified EIF5A2 protein (Cat# TP306249). The protein was produced from HEK293T cells transfected with EIF5A2 cDNA clone (Cat# [RC206249]) using MegaTran 2.0 (Cat# [TT210002]).