

Product datasheet for TP306269

SCHIP1 (NM_014575) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human schwannomin interacting protein 1 (SCHIP1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206269 protein sequence Red =Cloning site Green =Tags(s)
	MERSGQRVTTWDCDQGKHSDSDYREDGMDLGSDAGSSSSSRASSQSNSTKVTPCSECKSSSSPGGSLDL VSALEDYEEPFPVYQKKVIDEWAPEEDGEEEEEDERDQRGYRDDRSPAREPGDVSARTRSGGGGGRSAT TAMPPPVPNGNLHQHDPQDLRHNGNWWVAGRPCSRGRRAIKQPQAGRRSRGPAAGGLCLQPPDGGT CVPEEPPVPPMDWEALEKHLAQLFREQEVNRNQGQARTNSTSAQKNERESIRQKLALGFFDDGPGIYTS CSKSGKPSLSSRLQSGMNLQICFVNDSGSDKSDADDSKTETSLDTPLSPMSKQSSSYSDRDTTEEESES LDDMDFLTRQKKLQAEAKMALAMAKPMAKMQVEVEKQNRKKSPVADLLPHMPHISECLMKRSLKPTDLRD MTIGQLQVIVNDLHSQIESLNEELVQLLLIRDELHTEQDAMLVDIEDLTRHAESQQKHMAEKMPAK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	53.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_055390



[View online »](#)

Locus ID: 29970
UniProt ID: [Q9P0W5](#), [P0DPB3](#)
RefSeq Size: 2647
Cytogenetics: 3q25.32-q25.33
RefSeq ORF: 1392
Synonyms: SCHIP-1

Product images:



Coomassie blue staining of purified SCHIP1 protein (Cat# TP306269). The protein was produced from HEK293T cells transfected with SCHIP1 cDNA clone (Cat# [RC206269]) using MegaTran 2.0 (Cat# [TT210002]).