

Product datasheet for TP306300

GABPA (NM_002040) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human GA binding protein transcription factor, alpha subunit 60kDa (GABPA), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206300 protein sequence Red =Cloning site Green =Tags(s)

MTKREAELIEIEIDGTEKAECTEESIVEQTYAPAECVSQAIDINEPIGNLKKLLEPRLQCSLDAHEICL
QDIQLDPERSLFDQGVKTDGTVQLSVQVISYQGIEPKLNILEIVKPADTVEVIDPPDAHHAESEAHLVEE
AQVITLDGTKHITTISDETSEQVTRWAAALEGYRKEQERLGIPYDPIQWSTDQVLHWWVWVMKEFSMTDI
DLTTLNISGRELCSLNQEDFFQRVPRGEILWSHLELLRKYVLASQEQMNEIVTIDQPVQIIPASVQSAT
PTTIKVINSSAKAAKVQRAPRISGEDRSSPGNRTGNNGQIQLWQFLLELLTDKDARDDCISWVGDEGEFKL
NQPELVAQKWGQRKNKPTMNYEKLRSALRYYYDGDGMICKVQGKRFFVYKFVCDLKTLLIGYSAAELNRLVTE
CEQKKLAKMQLHGIAQPVTAVALSTASLQTEKDN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	51.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_002031](#)

Locus ID: 2551

UniProt ID: [Q06546](#), [A8IE48](#), [Q8IYS3](#)

RefSeq Size: 5182

Cytogenetics: 21q21.3


RefSeq ORF: 1362

Synonyms: E4TF1-60; E4TF1A; NFT2; NRF2; NRF2A; RCH04A07

Summary: This gene encodes one of three GA-binding protein transcription factor subunits which functions as a DNA-binding subunit. Since this subunit shares identity with a subunit encoding the nuclear respiratory factor 2 gene, it is likely involved in activation of cytochrome oxidase expression and nuclear control of mitochondrial function. This subunit also shares identity with a subunit constituting the transcription factor E4TF1, responsible for expression of the adenovirus E4 gene. Because of its chromosomal localization and ability to form heterodimers with other polypeptides, this gene may play a role in the Down Syndrome phenotype. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Oct 2010]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified GABPA protein (Cat# TP306300). The protein was produced from HEK293T cells transfected with GABPA cDNA clone (Cat# [RC206300]) using MegaTran 2.0 (Cat# [TT210002]).