

Product datasheet for TP306313

RALYL (NM_173848) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RALY RNA binding protein-like (RALYL), transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206313 protein sequence Red=Cloning site Green=Tags(s)

MTGKTQTSNVTNKNPKSINSRVFIGNLNTAIVKKVDIEAIFSKYGKIVGCSVHKGYAFVQYMSERHARA
AVAGENARVIAGQPLDINMAGEPKPYRPKPGNKRPLSALYRLESKEPFLSVGGYVFDYDYRDDFYNRFLF
DYHGRVPPPPRAVIPLKRPRVAVTTRRGKGVFSMKGGSRSTASGSTGSKLKSDELQTIKKELTQIKTKI
DSLLGRLEKIEKQQAEEAQKKQLEESLVLIQEECVSEIADHSTEEPAEGGPDADGEEMTDGIEEDFDE
DGGHELFLQIK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	32.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_776247
Locus ID:	138046



[View online »](#)

UniProt ID:	Q86SE5
RefSeq Size:	2919
Cytogenetics:	8q21.2
RefSeq ORF:	873
Synonyms:	HNRPCL3

Product images:

Coomassie blue staining of purified RALYL protein (Cat# TP306313). The protein was produced from HEK293T cells transfected with RALYL cDNA clone (Cat# [RC206313]) using MegaTran 2.0 (Cat# [TT210002]).