

Product datasheet for **TP306317L**

RPL18 (NM_000979) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human ribosomal protein L18 (RPL18), 1 mg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC206317 protein sequence
Red=Cloning site **Green**=Tags(s)

MGVDIRHNKDRKVRKPKSQDIYLRLLVKLYRFLARRTNSTFNQVVLKRLFMSRTNRPPLCLSRMIRKM
KLPGRENKTAVVGTITDDVRVQEVPKLKVCALRVTSRARSRILRAGGKILTFDQLALDSPKGCCTVLLS
GPRKGREVYRHF GKAPGTPHSHTKPYVRSKGRKFERARGRRASRGYKN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 21.5 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_000970](#)

Locus ID: 6141

UniProt ID: [Q07020](#), [A0A024QZD1](#)

RefSeq Size: 893



[View online »](#)

Cytogenetics: 19q13.33

RefSeq ORF: 564

Synonyms: DBA18; L18

Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a member of the L18E family of ribosomal proteins that is a component of the 60S subunit. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]

Protein Pathways: Ribosome

Product images:



Coomassie blue staining of purified RPL18 protein (Cat# [TP306317]). The protein was produced from HEK293T cells transfected with RPL18 cDNA clone (Cat# [RC206317]) using MegaTran 2.0 (Cat# [TT210002]).