

Product datasheet for **TP306350L**

SUPV3L1 (NM_003171) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human suppressor of var1, 3-like 1 (<i>S. cerevisiae</i>) (SUPV3L1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206350 protein sequence Red =Cloning site Green =Tags(s)

MSFSRALLWARLPAGRQAGHRAAICSA LRPHFGPFPGVLGQVSVLATASSASGGSKIPNTSLFVPLTVK
PQGPSADGDVGAE LTRPLDKNEVKKVLDFYKRKEIQKL GADYGLDARLFHQAFISFRNYIMQSHSLDVD
IHIVLNDICFGAAHADDLPFFLRHAKQIFPVLDCDDLRKISDLRIPPNWYPDARAMQRKIIFHSGPTN
SGKTYHAIQKYFSAKSGVYCGPLKLLAHEIFEKSNAAGVPCDLVTGEERTVQPNGKQASHVSCTVEMCS
VTPYEVAVIDEIQMIRDPARGWAWTRALLGLCAEEVHLCGEPA AIDLVMELMYTTGEEVEVRDYKRLTP
ISVLDHALES LDNL R PGDCIVCF SKNDIYSVSRQIEIRGLES AVIYGS LPPG TKLAQAKKFNDPNDPCKI
LVATDAIGMGLNLSIRRIIFYSLIKPSINEKGERELEPITTSQALQIAGRAGRFSSRFKEGEVTTMNHED
LSLLKEILKRPVDPIRAAGLHPTAEQIEMFAYHLPDATLSNLIDIFVDFSQVDGQYFVCNMDDFKFSAEL
IQHIPLSLRVRYVCTAPINKKQPFVCSSLLQFARQYSRNEPLTFAWLRRYIKWPLLPKNIKDLMDLEA
VHDVLDLYLWLSYRFMDMFDPASLIRDLQKELDGIIQDGVHNITKLIK MSETHKLLNLEGFPSGSQSRLS
GTLKSQARRTRGTKALGSKATEPPSPDAGELSLASRLVQQGLLTPDMLKQLEKEWMTQQTEHNKEKTESG
THPKGTRRRKKKEPDSD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	87.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003162
Locus ID:	6832
UniProt ID:	Q8IYB8
RefSeq Size:	2530
Cytogenetics:	10q22.1
RefSeq ORF:	2358
Synonyms:	SUV3
Summary:	Major helicase player in mitochondrial RNA metabolism. Component of the mitochondrial degradosome (mtEXO) complex, that degrades 3' overhang double-stranded RNA with a 3'-to-5' directionality in an ATP-dependent manner. ATPase and ATP-dependent multisubstrate helicase, able to unwind double-stranded (ds) DNA and RNA, and RNA/DNA heteroduplexes in the 5'-to-3' direction. Plays a role in the RNA surveillance system in mitochondria; regulates the stability of mature mRNAs, the removal of aberrantly formed mRNAs and the rapid degradation of non coding processing intermediates. Also implicated in recombination and chromatin maintenance pathways. May protect cells from apoptosis. Associates with mitochondrial DNA.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified SUPV3L1 protein (Cat# [TP306350]). The protein was produced from HEK293T cells transfected with SUPV3L1 cDNA clone (Cat# [RC206350]) using MegaTran 2.0 (Cat# [TT210002]).