

Product datasheet for **TP306358**

TTC9B (NM_152479) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tetratricopeptide repeat domain 9B (TTC9B), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206358 protein sequence Red =Cloning site Green =Tags(s)
	 MQRGALSPVLMLSAAPEPPRPPPALSPPGSGPGSGSRHGSARPGTPEPSGSLGAALDSSLRAAVAFKA EGQRCYREKKFREAIQKYHRALLQLKAAQGARPSGLPAPAPGPTSSPGPARLSEEQRRRLVESTEVECYDS LTACLLQSELVNYERVREYCLKVLEKQQGNFKATYRAGIAFYHLGDYARALRYLQEARREPTDTNVLRY IQLTQLKMNRCSLQREDSGAGSQTRDVIG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	25.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_689692</u>
Locus ID:	148014
UniProt ID:	<u>Q8N6N2</u>




[View online »](#)

RefSeq Size: 859

Cytogenetics: 19q13.2

RefSeq ORF: 717

Product images:



Coomassie blue staining of purified TTC9B protein (Cat# TP306358). The protein was produced from HEK293T cells transfected with TTC9B cDNA clone (Cat# [RC206358]) using MegaTran 2.0 (Cat# [TT210002]).