

Product datasheet for **TP306358M**

TTC9B (NM_152479) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human tetratricopeptide repeat domain 9B (TTC9B), 100 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC206358 protein sequence
Red=Cloning site **Green**=Tags(s)

MQRGALSPVLMLSAAPEPPRPPPALSPPGSGPGSGSRHGSARPGTPEPSGSLGAALDSSLRAAVAFKA
EGQRCYREKKFREAIQKYHRALLQLKAAQGARPSGLPAPAPGPTSSPGPARLSEEQRRRLVESTEVECYDS
LTACLLQSELVNYERVREYCLKVLEKQQGNFKATYRAGIAFYHLGDYARALRYLQARSREPTDTNVLRY
IQLTQLKMNRCSLQREDSGAGSQTRDVIG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 25.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_689692](#)

Locus ID: 148014

UniProt ID: [Q8N6N2](#)




[View online »](#)

RefSeq Size: 859

Cytogenetics: 19q13.2

RefSeq ORF: 717

Product images:



Coomassie blue staining of purified TTC9B protein (Cat# [TP306358]). The protein was produced from HEK293T cells transfected with TTC9B cDNA clone (Cat# [RC206358]) using MegaTran 2.0 (Cat# [TT210002]).