

Product datasheet for **TP306387L**

FAM71F1 (NM_032599) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 71, member F1 (FAM71F1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206387 protein sequence Red =Cloning site Green =Tags(s)

MLSSFPHRKTWRKSKKTVKVTRSYPTFPSLNAWEEFRGLLPVDGEPNPGAGLGVEEGLLCRVVHSPEFNL
FLDSVVFESNFIQVKRGRNWRDVYKASNTMALGVTSSVPCPLPNILLMASVKWHQGQNQTWNRPSIAPN
IFLKRIPLRFVELQVCDHYQRILQLRTVTEKIYYLKLHPDHPETVHFHWIRLVQILQKGLSITTKDPRI
LVTHCLVPKNCSSPSGDSKLVQKKLQASQPSESLIQLMTKGESEALSQIFADLHQQNQLRSSRKVETNKN
SSGKDSREDSIPTCDLRWRASFTYGWERENPSGLQPLSLLSTLAASTGPQLAPPIGNSI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	38.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_115988
Locus ID:	84691



[View online »](#)

UniProt ID: [Q96KD3](#), [A0A140VJJ3](#)
RefSeq Size: 1781
Cytogenetics: 7q32.1
RefSeq ORF: 3045
Synonyms: FAM137A; NYD-SP18

Product images:



Coomassie blue staining of purified FAM71F1 protein (Cat# [TP306387]). The protein was produced from HEK293T cells transfected with FAM71F1 cDNA clone (Cat# [RC206387]) using MegaTran 2.0 (Cat# [TT210002]).